



東京大学
THE UNIVERSITY OF TOKYO

Development of an absorption bench and some results about crystalline coatings

Manuel Marchiò, Raffaele Flaminio,
Danièle Forest, Laurent Pinard,
Daisuke Tatsumi, Kiichi Sakai




Kagra Face2Face Meeting
2017 年 3 月 30 日




Motivation

- Study the absorption of Sapphire, scans and maps
- Characterize the KAGRA mirrors, bulk and surface absorption
- R&D – Investigate on Crystalline coating to reduce mechanical losses.



Absorption measurement method at NAOJ




- The 1064nm pump periodically changes the temperature of the sample
- The thermal lens effect changes the probe wavefront
- After some propagation, the perturbation makes interference with the main mode
- The detector sense the intensity variation of the central part of the spot.
- This variation is proportional to the pump **absorption rate**

Photo-thermal
Common-path
Interferometer

from the company
Stanford Photo-Thermal Solution




Experimental setup



- Scan along the sample depth
- 2D maps of the surface and inside the substrate
- High power pump laser (25W)
- Best sensitivity:
 - Silica: <1 ppm/cm (bulk) <0.1 ppm (surface)
 - Sapphire: <5 ppm/cm (bulk) <0.5 ppm (surface)

Experimental setup




- Scan along the sample depth
- 2D maps of the surface and inside the substrate
- High power pump laser (25W)
- Best sensitivity:
 - Silica: <1 ppm/cm (bulk) <0.1 ppm (surface)
 - Sapphire: <5 ppm/cm (bulk) <0.5 ppm (surface)


Experimental setup

Large Translation Stage


KAGRA mirror size



TAMA mirror size



Up to 2" mirror size






Experimental setup

LABVIEW SOFTWARE

Commercial software couldn't control the large translation stage


→ We made it from scratch



Experimental setup




- Clean booth in TAMA300 central room
- HEPA SS-MAC-103 (2 units)
- Currently:
~7 particles/L (0.3 μ m)



Want it cleaner!


- Will replace the filter with new one
- Will remove all unnecessary objects and clean up often.
- Will use clean suits and prepare a clean entrance.



Calibration

- Reference sample used now:
Schott glass NG12
3.6mm thick
abs=22%/cm (spectrophotometer)

- For thicker samples
→ Photo-detector position correction
 $\text{Thickness} * (n-1) / n$





from presentation JPS autumn 2016
JGW-G1605695-v1 for more details

- Different materials have different thermal response
→ calibration correction factor is needed
NG12 → silica **0.73**
NG12 → sapphire **3.34**
provided by the company Stanford Photo-Thermal Solution

Calibration


We measured TAMA300 size sapphire samples



Brought it to LMA in Lyon (France)
use Photodeflection technique
Measured value = **25 ppm/cm**

We have to improve the
calibration accuracy

PHOTOTHERMAL DEFLECTION METHOD



Calibration


- To avoid material correction factors
→ will use a Titanium doped sapphire to calibrate
Low doping rate (0.01%) → same thermal proprieties
Absorption @1064 ~2 %/cm
We can measure it with spectrophotometer



Ti:Sapphire
Abs: ~2%/cm

- In addition, collaboration with Mio-san
use **Michelson absorption method** at Mio laboratory
for calibration check

R&D on Crystalline coatings



- Grown with Molecular Beam Epitaxy (MBE) on a **GaAs** substrate
- Then transferred onto **silica** or **sapphire** substrate


Optical performances

Measurements at LMA (France)	Coating on silica substrate	Coating on sapphire substrate
Transmission @1064nm	6ppm	6ppm
Roughness	1nm RMS	0.1nm RMS
Scattering @1064nm	9.5ppm	6ppm
Absorption @1064nm (1W)	0.4 ± 0.4 ppm/cm	<below noise>

Presentation JGW-G1605695-v1 for more details



R&D on Crystalline coatings



- Grown with Molecular Beam Epitaxy (MBE) on a **GaAs** substrate
- Then transferred onto **silica** or **sapphire** substrate

Optical performances


Measurements at LMA (France)	Coating on silica substrate	Coating on sapphire substrate
Transmission @1064nm	6ppm	6ppm
Roughness	1nm RMS	0.1nm RMS
Scattering @1064nm	9.5ppm	6ppm
Absorption @1064nm (1W)	0.4 ± 0.4 ppm/cm	<below noise>

Presentation JGW-G1605695-v1 for more details




R&D on Crystalline coatings

- Currently installing a 1310nm probe laser



- We will be able to measure GaAs samples, which are not transparent to 633nm
- We have a **25W power pump** laser (1064nm), more precision
- We have a translation stage to make maps of the surface absorption



Design of two probe lasers (633nm and 1310nm)



Summary


- Assembly of large translation stage ✓
- Able to measure thick sample ✓
- LabView software for automatic scan and maps ✓
- R&D: good results on crystalline coatings ✓
- Soon make the environment perfect clean
- Soon improve the calibration accuracy
- Soon be able to measure KAGRA sapphire test masses
- Soon be able to measure GaAs samples

Summary

- Assembly of large translation stage ✓
- Able to measure thick sample ✓
- LabView software for automatic scan and maps ✓
- R&D: good results on crystalline coatings ✓
- Soon make the environment perfect clean
- Soon improve the calibration accuracy
- Soon be able to measure KAGRA sapphire test masses
- Soon be able to measure GaAs samples

Thank you for your attention

The probe angle in the material (Snell's law)



$$\sin(\alpha_0)n_0 = \sin(\alpha_1)n_1$$

$$X_1 = X_0 \frac{\tan(\alpha_0)}{\tan(\alpha_1)} = \frac{n_1 \cos(\alpha_1)}{n_0 \cos(\alpha_0)}$$

$$n_0 = 1$$

$$n_1 = 1.76 \text{ Sapphire}$$

$$\alpha_0 = 6^\circ$$

$$X_1 = X_0 \cdot 1.81$$

To make a scan surface-to-surface of a KAGRA substrate we have to translate the sample by $150\text{mm}/1.81 = 83 \text{ mm}$ (instead of 150mm)

This is important in the design the new setup